

MOGAT1 Antibody



PACO10543

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MOGAT1

Uniprot

Q96PD6

Synonyms:

monoacylglycerol O-acyltransferase 1; MOGAT1; DGAT2L; DGAT2L1; MGAT1 ;

Immunogen:

Human MOGAT1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A