

# MPHOSPH8 Antibody



PACO10553

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

MPHOSPH8

**Uniprot**

Q99549

**Synonyms:**

M-phase phosphoprotein 8; MPHOSPH8; FLJ35237; HSMPP8; TWA3; mpp8 ;

**Immunogen:**

Human MPHOSPH8.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A