

# MPP6 Antibody



PACO10557

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MPP6

**Uniprot**

Q9NZW5

**Synonyms:**

membrane protein, palmitoylated 6 (MAGUK p55 subfamily member 6); MPP6; PALS2; VAM-1; VAM1; p55T ;

**Immunogen:**

Human MPP6.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A