

MPV17 Antibody



PACO10562

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MPV17

Uniprot

P39210

Synonyms:

MpV17 mitochondrial inner membrane protein; MPV17; SYM1 ;

Immunogen:

Human MPV17.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A