

MPZL1 Antibody



PACO10563

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MPZL1

Uniprot

O95297

Synonyms:

myelin protein zero-like 1; MPZL1; FLJ21047; PZR; PZR1b; PZRa; PZRb ;

Immunogen:

Human MPZL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A