

MRPL10 Antibody



PACO10572

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MRPL10

Uniprot

Q7Z7H8

Synonyms:

mitochondrial ribosomal protein L10; MRPL10; L10MT; MGC17973; MRP-L10; MRP-L8; MRPL8; RPML8 ;

Immunogen:

Human MRPL10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A