

MRPL38 Antibody



PACO10588

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

MRPL38

Uniprot

Q96DV4

Synonyms:

mitochondrial ribosomal protein L38; MRPL38; FLJ13996; HSPC262; KIAA1863; MGC4810; MRP-L3; RPML3 ;

Immunogen:

Human MRPL38.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A