

MRPL9 Antibody



PACO10595

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

MRPL9

Uniprot

Q9BYD2

Synonyms:

mitochondrial ribosomal protein L9; MRPL9; L9mt ;

Immunogen:

Human MRPL9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A