

# MRPS16 Antibody



PACO10600

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MRPS16

**Uniprot**

Q9Y3D3

**Synonyms:**

mitochondrial ribosomal protein S16; MRPS16; CGI-132; COXPD2; FLJ22062; FLJ40972; MRP-S16; RPMS16 ;

**Immunogen:**

Human MRPS16.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A