

# MRPS18A Antibody



PACO10602

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MRPS18A

**Uniprot**

Q9NVS2

**Synonyms:**

mitochondrial ribosomal protein S18A; MRPS18A; FLJ10548; HumanS18b; MRP-S18-3; MRPS18-3; S18bmt ;

**Immunogen:**

Human MRPS18A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A