

MRPS22 Antibody



PACO10606

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MRPS22

Uniprot

P82650

Synonyms:

mitochondrial ribosomal protein S22; MRPS22; C3orf5; COXPD5; GIBT; GK002; MRP-S22; RPMS22 ;

Immunogen:

Human MRPS22.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A