

MRPS24 Antibody



PACO10608

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MRPS24

Uniprot

Q96EL2

Synonyms:

mitochondrial ribosomal protein S24; MRPS24; HSPC335; MRP-S24 ;

Immunogen:

Human MRPS24.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A