

MRPS26 Antibody



PACO10610

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MRPS26

Uniprot

Q9BYN8

Synonyms:

mitochondrial ribosomal protein S26; MRPS26; C20orf193; GI008; MRP-S13; MRP-S26; MRPS13; NY-BR-87; RPMS13; dJ534B8.3 ;

Immunogen:

Human MRPS26.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A