

# MRPS27 Antibody



PACO10611

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MRPS27

**Uniprot**

Q92552

**Synonyms:**

mitochondrial ribosomal protein S27; MRPS27; FLJ21764; FLJ23348; KIAA0264; MRP-S27; S27mt ;

**Immunogen:**

Human MRPS27.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A