

MRPS34 Antibody



PACO10614

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MRPS34

Uniprot

P82930

Synonyms:

mitochondrial ribosomal protein S34; MRPS34; MGC2616; MRP-S12; MRP-S34; MRPS12 ;

Immunogen:

Human MRPS34.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A