

MS4A12 Antibody



PACO10619

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MS4A12

Uniprot

Q9NXJ0

Synonyms:

membrane-spanning 4-domains, subfamily A, member 12; MS4A12; Ms4a10 ;

Immunogen:

Human MS4A12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A