

# MSH2 Antibody



PACO10622

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MSH2

**Uniprot**

P43246

**Synonyms:**

mutS homolog 2, colon cancer, nonpolyposis type 1 (E. coli); MSH2; COCA1; FCC1; HNPCC; HNPCC1; LCFS2 ;

**Immunogen:**

Human MSH2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A