

MSL3 Antibody



PACO10624

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MSL3

Uniprot

Q8N5Y2

Synonyms:

male-specific lethal 3 homolog (Drosophila); MSL3; DKFZp586J1822; MSL3L1 ;

Immunogen:

Human MSL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A