

MSRB2 Antibody



PACO10628

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MSRB2

Uniprot

Q9Y3D2

Synonyms:

methionine sulfoxide reductase B2; MSRB2; CBS-1; CBS1; CGI-131; MGC26104; MSRB; PILB ;

Immunogen:

Human MSRB2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A