

MTHFD1L Antibody



PACO10639

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

MTHFD1L

Uniprot

Q6UB35

Synonyms:

methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 1-like; MTHFD1L; DKFZp586G1517; FLJ21145; FTHFSDC1; MTC1THFS; dJ292B18.2 ;

Immunogen:

Human MTHFD1L

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A