

MTIF3 Antibody



PACO10642

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MTIF3

Uniprot

Q9H2K0

Synonyms:

mitochondrial translational initiation factor 3; MTIF3; IF-3mt; IF3(mt) ;

Immunogen:

Human MTIF3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A