

MTMR12 Antibody



PACO10644

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MTMR12

Uniprot

Q9C0I1

Synonyms:

myotubularin related protein 12; MTMR12;3-PAP; PIP3AP ;

Immunogen:

Human MTMR12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A