

# FAN1 Antibody



PACO10646

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FAN1

**Uniprot**

Q9Y2M0

**Synonyms:**

myotubularin related protein 15; MTMR15; DKFZp451H236; DKFZp686K16147; KIAA1018;

**Immunogen:**

Human MTMR15.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A