

MTMR4 Antibody



PACO10648

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MTMR4

Uniprot

Q9NYA4

Synonyms:

myotubularin related protein 4; MTMR4; FYVE-DSP2; KIAA0647; ZFYVE11 ;

Immunogen:

Human MTMR4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A