

MTX2 Antibody



PACO10656

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MTX2

Uniprot

O75431

Synonyms:

metaxin 2; MTX2; MGC111067 ;

Immunogen:

Human MTX2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A