

MUM1 Antibody



PACO10659

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MUM1

Uniprot

Q2TAK8

Synonyms:

melanoma associated antigen (mutated) 1; MUM1; FLJ14868; FLJ22283; HSPC211; MGC131891; MGC163315; MUM-1 ;

Immunogen:

Human MUM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A