

MYB Antibody



PACO10668

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYB

Uniprot

P10242

Synonyms:

v-myb myeloblastosis viral oncogene homolog (avian); MYB; Cmyb; c-myb; c-myb_CDS; efg ;

Immunogen:

Human MYB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A