

MYBBP1A Antibody



PACO10669

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYBBP1A

Uniprot

Q9BQG0

Synonyms:

MYB binding protein (P160) 1a; MYBBP1A; FLJ37886; P160; PAP2 ;

Immunogen:

Human MYBBP1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A