

MYBL2 Antibody



PACO10670

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYBL2

Uniprot

P10244

Synonyms:

v-myb myeloblastosis viral oncogene homolog (avian)-like 2; MYBL2; B-MYB; BMYB; MGC15600 ;

Immunogen:

Human MYBL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A