

MYEF2 Antibody



PACO10676

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYEF2

Uniprot

Q9P2K5

Synonyms:

myelin expression factor 2; MYEF2; FLJ11213; HsT18564; KIAA1341; MEF-2; MGC87325; MST156; MSTP156 ;

Immunogen:

Human MYEF2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A