

MYF6 Antibody



PACO10678

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MYF6

Uniprot

P23409

Synonyms:

myogenic factor 6 (herculin); MYF6; MRF4; bHLHc4 ;

Immunogen:

Human MYF6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A