

MYH10 Antibody



PACO10679

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYH10

Uniprot

P35580

Synonyms:

myosin, heavy chain 10, non-muscle; MYH10; MGC134913; MGC134914; NMMHCB ;

Immunogen:

Human MYH10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A