

MYL1 Antibody



PACO10684

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MYL1

Uniprot

P05976

Synonyms:

myosin, light chain 1, alkali; skeletal, fast; MYL1; MLC1F; MLC3F ;

Immunogen:

Human MYL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A