

MYL2 Antibody



PACO10686

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MYL2

Uniprot

P10916

Synonyms:

myosin, light chain 2, regulatory, cardiac, slow; MYL2; CMH10; DKFZp779C0562; MLC2 ; MLC-2v; MLC-2

Immunogen:

Human MYL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A