

MYL4 Antibody



PACO10688

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MYL4

Uniprot

P12829

Synonyms:

myosin, light chain 4, alkali; atrial, embryonic; MYL4; ALC1; AMLC; GT1; PRO1957 ;

Immunogen:

Human MYL4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A