

MYL5 Antibody



PACO10689

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYL5

Uniprot

Q02045

Synonyms:

myosin, light chain 5, regulatory;

Immunogen:

Human MYL5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A