

MYL9 Antibody



PACO10691

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYL9

Uniprot

P24844

Synonyms:

myosin, light chain 9, regulatory; MYL9; LC20; MGC3505; MLC2; MRLC1; MYRL2 ;

Immunogen:

Human MYL9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A