

MYOF Antibody



PACO10709

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYOF

Uniprot

Q9NZM1

Synonyms:

myoferlin; MYOF; FER1L3; FLJ36571; FLJ90777 ;

Immunogen:

Human MYOF.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A