

# MYOM3 Antibody



PACO10710

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MYOM3

**Uniprot**

Q5VTT5

**Synonyms:**

myomesin family, member 3; MYOM3; FLJ35961; MGC71483; RP11-293P20.1 ;

**Immunogen:**

Human MYOM3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A