

# NAMPT Antibody



PACO10722

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NAMPT

**Uniprot**

P43490

**Synonyms:**

nicotinamide phosphoribosyltransferase; NAMPT;1110035O14Rik; DKFZp666B131; MGC117256; PBEF; PBEF1; VF; VISFATIN ;

**Immunogen:**

Human NAMPT.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A