

NAP1L1 Antibody



PACO10726

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NAP1L1

Uniprot

P55209

Synonyms:

nucleosome assembly protein 1-like 1; NAP1L1; FLJ16112; MGC23410; MGC8688; NAP1; NAP1L; NRP ;

Immunogen:

Human NAP1L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A