

PACO10736

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NASP

Uniprot

P49321

Synonyms:

nuclear autoantigenic sperm protein (histone-binding); NASP; DKFZp547F162; FLB7527; FLJ31599; FLJ35510; MGC19722; MGC20372; MGC2297; PRO1999 ;

Immunogen:

Human NASP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A