

NAV1 Antibody



PACO10743

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NAV1

Uniprot

Q8NEY1

Synonyms:

neuron navigator 1; NAV1; DKFZp781D0314; FLJ12560; FLJ14203; KIAA1151; MGC14961; POMFIL3; steerin-1 ;

Immunogen:

Human NAV1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A