

# NBL1 Antibody



PACO10746

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NBL1

**Uniprot**

P41271

**Synonyms:**

neuroblastoma, suppression of tumorigenicity 1; NBL1; D1S1733E; DAN; DAND1; MGC8972; NB; NO3 ;

**Immunogen:**

Human NBL1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A