

PACO10747

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NBN

**Uniprot**

O60934

**Synonyms:**

nibrin; NBN; AT-V1; AT-V2; ATV; FLJ10155; MGC87362; NBS; NBS1; P95 ;

**Immunogen:**

Human NBN.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A