

NBR1 Antibody



PACO10748

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

NBR1

Uniprot

Q14596

Synonyms:

neighbor of BRCA1 gene 1; NBR1;1A1-3B; FLJ55359; FLJ98272; KIAA0049; M17S2; MIG19 ;

Immunogen:

Human NBR1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A