

NDE1 Antibody



PACO10768

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NDE1

Uniprot

Q9NXR1

Synonyms:

nudE nuclear distribution gene E homolog 1 (A. nidulans); NDE1; FLJ20101; HOM-TES-87; NUDE; NUDE1 ;

Immunogen:

Human NDE1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A