

NDUFAF1 Antibody



PACO10782

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NDUFAF1

Uniprot

Q9Y375

Synonyms:

NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, assembly factor 1; NDUFAF1; CGI-65; CGI65; CIA30 ;

Immunogen:

Human NDUFAF1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A