

NDUFB3 Antibody



PACO10788

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NDUFB3

Uniprot

O43676

Synonyms:

NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 3, 12kDa; NDUFB3; B12 ;

Immunogen:

Human NDUFB3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A