

NDUFV2 Antibody



PACO10800

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NDUFV2

Uniprot

P19404

Synonyms:

NADH dehydrogenase (ubiquinone) flavoprotein 2, 24kDa;

Immunogen:

Human NDUFV2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A